



SECOND SKIN PTY LTD 40 O'MALLEY STREET, OSBORNE PARK WA 6017

Existing Patient

New Patient

P: +61 8 9201 9455 E: orders@secondskin.com.au

or upload via www.secondskin.com.au/contact/enquiry (choose "an order")

PATIENT DETAILS FORM

Date:	New Order (✓)		Reorder (✓)									
PATIENT: (Surname)		(Given Names										
Date of Birth:					M 🗆	F□						
Patient Address:												
		Post Code:										
Patient Phone No: (Home)		(Work)										
HOSPITAL:		Order Number:										
Hospital Address:												
	Post Code:											
Therapist Name:		Department:										
Therapist Phone No:		Pager No:										
Therapist Email												
Photo Sent (✓) YES	NO	Email		POST/COL	JRIER							

GARMENT/GARMENTS REQUIRED:	
SEND ACCOUNT TO: (Include Claim/Reference Number)	
SEND GARMENT TO: Therapist - address as above (</td <td>Patient - address as above (✓)</td>	Patient - address as above (✓)
DATE REQUIRED BY:	

Second Skin will always endeavour to supply this order by the date you require.


Please keep in mind that delivery is subject to freight times and the receipt of written funding approval / hospital order numbers.

SECOND SKIN PTY LTD 40 O'Malley Street		E: orders@secondskin.com.au PAG	PAGE NO:			
OSBORNE PARK WA 6017						
CLIENT SURNAME: GI	VEN N	AME:	F 🗌 M DATE:	/	/	
Diagnosis: Burns 🗌 Lymphoedema 🗌 Trauma	ΠV	ascula	r Insufficiency 🗌 Other:			
Colour: Light Dark Black (Powersoft a						
Garment personalisation *Please choose carefully as garmer	nts canno	ot be exc	hanged/returned for change of mind or incorrect choice			
Stitching colour: (Circle one only) Purple/Green/Pink/Blu	e/Yello	w/Whi	te/Red/Orange			
Trim Colour: (Circle one only) Pink/Yellow/Green/Purple/	Navy/F	Red/Bla	ack/White			
Motif: (choose one only) M	otif co	lour: (choose one only)			
1. Style			7. Ankle	L	R	
Single leg			Centre front vertical seam (preferred option)			
Two leg			Ankle crease seam			
One and a half leg			Dorsal ankle gusset:			
Stump support			- Shimmer			
Panty girdle			- Powernet			
Flap tight			- Powersoft			
Hernia support			- Single hydrophobic			
Scrotal support			- Double hydrophobic			
All in one (see all in one form)			- Hydrophobic lining			
2. Fabric			8. Toes	L	R	
Powernet			Closed			
Powersoft			Big toe separate			
Shimmer			Foot glove			
Single hydrophobic			Stirrups			
Double hydrophobic			9. Zips - Lower Body	L	R	
3. Crotch			None in legs			
Open			Waist to thigh high			
Closed			Full length curved into foot			
Fly front			Below knee - straight medial to ankle			
4. Leg Lengths	L	R	Below knee - straight lateral to ankle			
Above knee			Below knee - curved medial into foot			
Ankle length			Below knee - curved lateral into foot		R	
Including feet 5a. Knee Gusset	L	R	10. Reinforcing Shimmer	L	ĸ	
Posterior knee gusset - shimmer	L	ĸ	Powernet			
Knee flexion gusset - all shimmer			Powersoft			
Knee flexion gusset - an similar			Sole			
Knee flexion gusset - powersoft anterior			Sole leather			
Knee flexion gusset - all single hydrophobic			Heel			
Knee flexion gusset - all double hydrophobic			Dorsum of foot			
5b. Hydrophobic Lining - Knee	I		Lower leg - anterior	+		
(a) anterior			Lower leg - posterior	-		
(b) posterior			Full leg - anterior	-		
(c) circumferential			Full leg - posterior	-		
6. Dressing Assist			11. Additional Options		1	
Zip tab			Colostomy site with hole and zip access			
Zip looper			Shaped abdomen			
Leather assist			Pregnancy panel			
			Soft braces with velcro closure			

Note any further design options you require. Call our design department in Perth (08 9201 9455) for any queries



PAGE NO: _____





PAGE NO:

FOOT TRACE FORM

CONFIDENTIAL

CLIE	ENT SURNAME	:				GIVE		≣:				. 🗌 F	🗌 M	DA	TE:	_//	/
 									0								
cm x 25 cm																	
Scale 1:1 19								LO ^K)						
Grid to							CONCE ON CE										
											Γ	1025	uripo	ı Tin	Q		
	Important: Measure leng trace from tip	ace from tip of big toe to tip of heel.							 Measuring Tips For big toe separate, measure big toe circumference and length. For a Foot Glove measure all toe circumferences and lengths Circumference measurements are taken at the middle of toe. Length measurements are taken from web space to tip of toe on the side of the toe as indicated with length arrow. 								