



Size and colour chart

Tights – Above Ankle – Child Open or Closed Pubis

COLOURS



FABRIC OPTIONS




Size	Waist (a)	Hips (b)	Top leg (c)	Knee (d)	Inside leg length to ankle (e)
2	47 – 50cm	47.5 – 51.5cm	26.5 – 31cm	20 – 21.5cm	34cm
4	50 – 53.5cm	51.5 – 57cm	31 – 35cm	21.5 – 23.5cm	40cm
6	53.5 – 57cm	57 – 63cm	35 – 38cm	23.5 – 25.5cm	46cm
8	57 – 60cm	63 – 69cm	38 – 42cm	25.5 – 27.5cm	52cm
10	60 – 64cm	69 – 76cm	42 – 46cm	27.5 – 29.5cm	58cm
12	64 – 68cm	76 – 83cm	46 – 50cm	29.5 – 32.5cm	64cm

MEASUREMENTS

Measuring Tips:

- All circumferences to be measured to the nearest millimetre.
- Tape should be snug NOT tight.
- Waist measurement is taken over the navel.
- Hip measurement is taken over top of hip bone.
- Knee measurement is over the middle of the patella.
- Leg length is measured to the end of the tights.



Please note that garments cannot be altered or repaired by Second Skin.



Second Skin Pty Ltd
40 O'Malley Street
Osborne Park WA Australia 6017
Email: perth@secondskin.com.au
Phone: +61 8 9201 9455

secondskin.com.au





Second Skin

CUSTOMISED COMPRESSION GARMENTS & SPLINTS

Fabric options

All Second Skin garments are made from a range of medical grade stretch fabrics. Several highly durable medical grade fabrics are incorporated in our garment designs, to address the individual needs of each client.

POWERNET FABRICS

Powernet is an open weave fabric which is prescribed for healed skin and for areas where additional support and compression is required. This fabric is particularly durable in work, school or outdoor activities. Powernet is often used in burns scar management and lymphoedema.

Our Powernet fabric now incorporates patented cooling properties, HeiQ Smart Temp, which provides outstanding thermal comfort to clients by keeping them cooler and dryer while wearing their garment.

Powernet fabric is available in light tan, dark tan and black.

SHIMMER FABRICS

Shimmer is a closed weave fabric which provides a softer interface with the skin. Shimmer fabric contours well and is comfortable and effective over joint surfaces and areas of skin fragility such as the elbow, ankle and knee. Shimmer fabric is used for early burns and flap management, external fixators, vascular conditions and certain types of lymphoedema.

Shimmer fabric is available in light tan, dark tan and black.

HYDROPHOBIC FABRICS

Hydrophobic fabric is used as both a lining and a base fabric in compression garments. This fabric wicks moisture away from the surface of the skin to evaporate, protecting the skin surface and ensuring optimal comfort.

Our hydrophobic fabric has an antimicrobial component embedded within the matrix of the fabric. This ensures our products are breathable and highly effective at removing bacteria from the fabric surface, to provide the wearer with all year round comfort. Due to this innovation, our garments are frequently prescribed to assist with the management of repeated skin breakdown and fragility.

Hydrophobic fabric is available in light tan, dark tan and black.

PERTH OFFICE

40 O'Malley Street
OSBORNE PARK WA 6017

P: +61 8 9201 9455

E: perth@secondskin.com.au

SYDNEY OFFICE

Suite 1, 8 Northcote Street
ST LEONARDS NSW 2065

P: +61 2 9386 0812

E: sydney@secondskin.com.au

BRISBANE OFFICE

Suite 1/41 Manilla Street
EAST BRISBANE QLD 4169

P: +61 7 3084 4319

E: brisbane@secondskin.com.au

MELBOURNE OFFICE

Unit 9, 202 Ferntree Gully Road
NOTTING HILL VIC 3168

P: +61 3 8637 0297

E: melbourne@secondskin.com.au